

# CellVision Counting Chambers

## Disposable cell counting slides

### Easy to use

**Fill-wait-count** disposable counting chambers with a clear and bright built-in grid for manual counting.

### Hygienic and efficient

You eliminate non-hygienic and time consuming cleaning procedures associated with standard counting chambers.


Your risk of being exposed to potentially infectious material is reduced drastically.

### Reliable counts


The CellVision slides are produced with a precise fixed depth with superior accuracy and reproducibility of the results.




### Built-in counting grids




10<sup>2</sup> grid




Neubauer improved



Bürker



Bürker-Türk



Thoma

Ref	Product	Package
<b>Counting chambers with grid</b>		
50001	CellVision 20 µm - 10 <sup>2</sup> grid (type Makler)	25 units/box
50002	CellVision 100 µm - Neubauer Improved grid	25 units/box
50003	CellVision 100 µm - Bürker grid	25 units/box
50004	CellVision 100 µm - Bürker-Türk grid	25 units/box
50005	CellVision 100 µm - Thoma grid	25 units/box
50006	CellVision 200 µm - Fuchs Rosenthal grid	25 units/box
<b>Counting chambers without grid</b>		
50105	CellVision 16 µm - 2 chambers (CASA)	25 units/box
50102	CellVision 20 µm - 2 chambers	25 units/box
50103	CellVision 20 µm - 4 chambers	25 units/box
50104	CellVision 100 µm - 2 chambers	25 units/box

