

Oosafe Disinfectants

For safely cleaning your entire IVF laboratory

The Green Line

Oosafe Incubator and LAF cabinets

An anticorrosive ready-to-use disinfectant for all laboratory equipment, including laminar flow hoods and the insides of incubators.



The Blue Line

Oosafe floor/laboratory surface

An economical (dilution 1:100) surface disinfectant containing a non-ionic soap. Dilute>apply>wait 15 min>rinse and wipe clean.




The Orange Line

Oosafe Hand


A ready-to-use hand disinfectant that will not dry out your skin.




Mouse embryos after 72 h incubation in disinfected MINC incubator




Control



6% hydrogen peroxide



70% ethanol



Oosafe®

Study done by:
William Lee, Emma Lingham,
Dr. Jim Catt, Dr. Sally Catt &
Dr. Mulyoto Pangestu

Department of Obstetrics
and Gynaecology

Education Centre of
Reproduction and
Development

MONASH University
Australia, 2012

- Microbial efficacy:
Against bacteria, fungus, and viruses
(hepatitis, HIV, influenza and rota virus)
- Active ingredient:
Quarternary ammonium compound
- No VOC release, no smell, no alcohol
- MEA-tested and CE-marked
- Long shelf-life:
Three years from manufacturing date,
even after opening

Ref	Product	Package
15550	Oosafe Incubator disinfectant	1 L
15551	Oosafe Incubator disinfectant	5 L
15557	Oosafe disinfectant wipes	70/bucket
15552	Oosafe Hand disinfectant	0,5 L
15553	Oosafe Hand disinfectant	5 L
15554	Oosafe Surface disinfectant	2 L
15555	Oosafe Surface disinfectant	10 L