

# TPC Micropipettes

Precision pipettes from 

## Superior quality



The TPC micropipettes are made by Australian embryologists who have been involved with reproductive medicine and micromanipulation techniques for many years.



Product	Tiplength	Od/Id	Ref 30°	Ref 35°
---------	-----------	-------	---------	---------

## Injection pipettes

Injection pipette w/spike	1.0 mm	7/5 µm	20100	20150
Injection pipette w/spike - extra rigid	1.0 mm	7/5 µm	20110	20160
Injection pipette w/spike - thin	0.6 mm	6.0/4.0 µm	20116	20166
Injection pipette w/spike - large (biopsies)	1.0 mm	8.0/6.0 µm	20113	20163
Injection pipette without spike 21°	0.6 mm	6.5/4.5 µm	20130	20180



## Holding pipettes

Holding pipette	1.0 mm	120/25 µm	20000	20050
Holding pipette - thin	1.0 mm	100/20 µm	20010	20060
Holding pipette - large	1.0 mm	120/30 µm	20012	20062



## Hatching pipettes

Hatching pipette	1.0 mm	12/10 µm	20300	20350
------------------	--------	----------	-------	-------



## Zona cutting pipettes (PZD)

Zona Cutting pipette	1.0 mm	2-3 µm	20400	20450
----------------------	--------	--------	-------	-------



## Biopsy pipettes

Biopsy pipette - bevel 40°	1.0 mm	30/28 µm	20200	20250
Biopsy pipette - straight 90°	1.0 mm	27/25 µm	20210	20260

