

UniHeat™ Table Plates

With integrated heating system

Table plates for IVF workstations

The table plates are custom made and can be of any size to replace existing table plates or to be incorporated into new cabinets.

Unique heating system

The liquid heating system is silent, free from vibrations and completely integrated into the table plate. Quality is kept at the highest level while the risk of contamination from open waterbaths is eliminated.

Complete solutions / separate parts

We supply complete IVF workstations or the individual parts like heated table plates, microscopes, gassing systems, cabinets and adjustable support stands, to help you keep expenditure to a minimum.

Ref	Products
70510	UniHeat™ table plate 120 cm - single
70511	UniHeat™ table plate 180 cm - single
70512	UniHeat™ table plate 180 cm - double
70513	UniHeat™ table plate - custom made

The standard UniHeat™ plates are designed for the classic Holten/K-Systems cabinet but can be made to fit any cabinet



The UniHeat™ heating system



The heating system is a completely closed liquid based heating system built into the table plate. It consists of two parts; an aluminum composite heat exchanger in the table plate, and a heater and a liquid circulator integrated into a liquid reservoir placed inside the table plate.

The heat uniformity is unique and by using a liquid based system we have eliminated cold spots and problems with overheating otherwise associated with electrically based heating systems.