

Brady BMP 51 label printer

New and improved design

Large display easy to operate

The extra-large display is very bright, clear and easy to read. The entire label is visible on the display making it easy to adjust the layout

Improved label cartridge

The design of the cartridge has been improved and issues with stuck labels and ripped ink ribbons no longer exist.


Large display easy to operate

The LabelMark Pro labeling software lets you easily design, edit and print labels. It is very intuitive, and serializing and multi-print are easily done.



Straw labels

Self laminating labels for cryo straws and vitrification straws.



FreezerBonz labels

Adhere perfectly to CryoTubes on clean, frozen or frosted surfaces.



Dish labels

White and transparent dish labels with printing under the dish.




Slide labels

Solvent resistant labels securing correct identification of the slides.



Ref	Product	Material	Type	A	B	C	QTY
72700	BMP51 label printer EU						1 printer
75502	LabelMark 5 Pro (software)	-	-	-	-	-	1 program
72710	Cryo straw label	M19427	I	25,4	25,4	9,53	200 labels
72711	Cryo straw label - short	M96427	I	25,4	19,05	9,53	200 labels
72712	Cryo label for CBS ID stick	M105427	III	27,6	6,35	-	360 labels
72720	White dish label	M81488	III	48,26	6,36	-	160 labels
72722	Slide label, solvent resistant	M97488	III	22,86	22,86	-	280 labels
72722	Transparent dish label	M17432	III	25,40	12,70	-	280 labels
72730	FreezerBonz label	M124490	III	41,91	12,7	-	180 labels



Note: Many more labels are available - specifications are sent on request