

# Craft™ Suction Pumps

## For oocyte aspiration



The Craft™ Suction Pump provides smooth low volume vacuum at pre-determined negative pressures, resulting in a low turbulence, oocyte recovery.

The duo-vac model includes a high vacuum function for clearing clutched aspiration needles etc.



### Single use fluid trap set

- replaces conventional reusable bottle set
- one patient one fluid trap set
- eliminates risk of contamination

Ref	Product
75000	Aspiration pump, CRAFT single
75010	Aspiration pump, CRAFT Duo-Vac
10412	Silicone tubing, 25 meters
10411	Male luerlock, polypropylene (25 units)
10410	Single Use Fluid Trap Set