

Cryo Cassettes

For convenient and safe storage of cryostraws

Less storage tanks are needed

The Cryocassettes allows more straws to be stored in each canister why the storage tanks are used more efficiently. Consequently less storage tanks, less nitrogen and less lab space is used.

Straws will never “float”

The straws are safely secured inside the cassettes making it impossible for the straws to float inside the storage tank.

Easy to use

Constructed with a convenient lifter making it easy and very fast to remove or insert a straw.



Capacity of the cryo cassettes

Diameter	Flat cassettes	Triangular cassettes
0.25 cc straws	8 straws / 3 with ID stick	24 straws / 12 with ID stick
0.5 cc straws	3 straws	12 straws

Wide range of colour coding options

The outer shell is transparent and the inner part comes in different colours for easy identification. Combined with colour coded straws you will have a very wide range of colour coding options making it easy to identify the individual straw.

Capacity of your canister

Diameter	Flat cassettes	Triangular cassettes
38 mm	12 units	6 units
56 mm	25 units	11 units
70 mm	37 units	18 units

Available for all canisters

The cryo cassettes are available in a flat version for storing embryos and a triangular version storing for sperm straws or bulk storage. Additionally, both versions come in a 140 mm and a 280 mm version to accommodate the canisters of your storage tank.



Color	Cassettes 140 mm		Cassettes 280 mm	
	Flat	Triangular	Flat	Triangular
Blue	57000	57020	57010	57030
Green	57001	57021	57011	57031
Red	57002	57022	57012	57032
Yellow	57003	57023	57013	57033
Beige	57004	57024	57014	57034
White	57005	57025	57016	57035