

Glass Pasteur Pipettes

MEA tested and CE marked for IVF use



IVF Glass Pasteur Pipettes

The pasteur pipettes are Mouse embryo tested and CE marked as a class IIa Medical Device under the EU medical directive.

Plugged and unplugged

The pipette comes in a short and a long version and a special synthetic material is used for the plugged version to avoid cotton wool fibre contamination.

Wide variety

The pipettes are made of either Borosilicate glass or soda glass, and packed in sleeves of 5 in boxes of 100 units.



Ref	Product	Package
30100	IVF Pasteur Pipette, 150 mm soda	100 units
30101	IVF Pasteur Pipette, 230 mm soda	100 units
30102	IVF Pasteur Pipette, 150 mm soda - plugged	100 units
30103	IVF Pasteur Pipette, 230 mm soda - plugged	100 units
30110	IVF Pasteur Pipette, 150 mm borosilicate	100 units
30111	IVF Pasteur Pipette, 230 mm borosilicate	100 units
30112	IVF Pasteur Pipette, 150 mm borosilicate - plugged	100 units
30113	IVF Pasteur Pipette, 230 mm borosilicate - plugged	100 units