

Soft HSU/HSG catheter

With separate stylet

Size
5 fr or 7 fr

For sonohysterography and hysterosalpingography

A soft latex-free balloon catheter with a rigid insertion sheath. The insertion sheath has a rounded tip and is a perfect guide and support during insertion eliminating the need for a tenaculum. The included syringe is preset with correct balloon volume - just inflate. For very difficult cases a stylet is available separately.

Ref	Product	Unit
13050	Zinnanti HS Catheter 5 fr	10/box
13051	Zinnanti HS Catheter 7 fr	10/box
13053	Shapeable HS Catheter 5 fr	10/box
13054	Shapeable HS Catheter 7 fr	10/box
13052	Stylet for Zinnanti HS Catheter	10/box

