

# IVF handling pipette

## Oocyte and embryo handling tip

**New**

### Correct and safe handling

Correct and safe handling of oocytes during pickup and when moving embryos from one culture media to another.

### Clear and non-breakable

The tip is made of PET which is crystal clear, non breakable and does not scratch dishes.

### CE marked and tested for IVF use

The tip is **CE** marked, **gamma** sterilized, **LAL** tested and **mouse embryo tested**.



Manufacturer: Gynetic Medical Products



Ref	Product	Packing
30012	Handling pipette 95 mm - single packed	50 units
30010	Handling pipette 55 mm - single packed	50 units
30011	Handling pipette 55 mm - 10 pack	250 units
15260	1 ml syringe - BD/Plastipack	100/box



**Tip Dimentions**  
 Total length: 55 & 95 mm  
 Inner diameter: 500 µm