

# Oosafe® Inline Filter

## For purifying incubator gas

### Reducing contaminants in incubator gas


Oosafe® filters are used inline between incubators and the gas sources in order to clean and purify the incoming gas.

### Protecting your gametes

Oosafe® filters protects your oocytes and embryos against harmful contaminants such as VOCs, CACs, toxins, bacteria, odor, dust and particles in the supplied incubator gas.

### Designed for IVF use

Oosafe® filters are CE marked and composed of specially designed active carbon filters and of high quality HEPA filters.



### Recommended usage

Oosafe® filters are single packed in boxes of 4 units and designed to trap contaminants for 90 days => recommended change every 3 months.

Ref	Product
73010	Oosafe® Inline Filter (4 units)
73011	Holder for Oosafe® filter
73012	Quick connector - Ø6 mm
73013	Quick connector - Ø8-10 mm
73014	Quick connector - Ø6 mm hard tubing
73015	Quick connector - Ø8 mm hard tubing

