

Speculum with side opening

Ideal for HSU/HSG examinations

Convenient side opening

The open side improves instrument access and allows the speculum to be easily removed over a catheter or other device that has been inserted into the uterus.

Left & Right designs

The speculum is available in 3 sizes with right side and left side opening, and it is made of durable quality German stainless steel.

Anatomically shaped

The lower blade is longer than the upper blade, making the cervix more easily visible.

Ref	Product
13080	Weisman-Graves Speculum - 7,5 x 2 cm, Left
13081	Weisman-Graves Speculum - 10 x 4 cm, Left
13082	Weisman-Graves Speculum - 12 x 4 cm, Left
13083	Weisman-Graves Speculum - 7,5 x 2 cm, Right
13084	Weisman-Graves Speculum - 10 x 4 cm, Right
13085	Weisman-Graves Speculum - 12 x 4 cm, Right