

Strip-ID pipettes

For embryo handling and oocyte denudation

Flexible and unbreakable


The Strip-ID pipette is a flexible unbreakable microcapillary used for handling oocytes and embryo during IVF procedures.

Clear and non-scratching

The Strip-ID pipette is crystal clear, sterile, MEA and LAL tested. It comes in a convenient 10-pack and is made of flexible polycarbonate which does not scratch culture dishes and embryo slides.

Universal fitting

The Strip-ID pipette is available in a variety of sizes and made with a compatible design to fit most commonly used pipette holders.



Ref	Product	Package
20550	Strip-ID pipette - 125	10/tube - 50 box
20551	Strip-ID pipette - 130	10/tube - 50 box
20552	Strip-ID pipette - 135	10/tube - 50 box
20553	Strip-ID pipette - 140	10/tube - 50 box
20554	Strip-ID pipette - 145	10/tube - 50 box
20555	Strip-ID pipette - 150	10/tube - 50 box
20556	Strip-ID pipette - 170	10/tube - 50 box
20557	Strip-ID pipette - 175	10/tube - 50 box
20558	Strip-ID pipette - 200	10/tube - 50 box

